MBH0063

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PRODUCT SHEET

Expired since 30/03/2023

IDENTIFICATION OF THE MEDICAL DEVICE	Type	Rechargeable battery	 
	Commercial designation	Batterie médicale rechargeable 10x V250H 10S1P ST5 12V 250mAh JST	
	Reference	MBH0063	
	EAN	3660766498175	
	Brand	NX	
	Compatible / Original battery	Compatible	
	Packaging	Unitary	

RECOMMENDED USAGE

Follow the instructions and recommendations specific to each model, using the technical instructions and document resources from the devices in which the battery is used

Brands	Equipment	Models	PN
Imex / Medasonics / Sonicaid	Medical device	CT Plus Doppler / Dop CT / Free Dop / D421 / 421 Fetal Doppler	5674/BR-002/BT0022/BT0024/C631

Identification	GENERAL TECHNICAL CHARACTERISTICS	Chemistry	NiMh
		Type	B250H
	IEC designation	10 HRM 6/25	
	Rated voltage	12V	
	Nominal capacity	250mAh	
	Internal resistance Ω	<500mΩ	

The voltage and the actual capacity in use can be affected by several factors, especially the temperature, the discharge current, the pack's history (ex:use, storage), etc

ELECTRICAL CHARACTERISTICS	CHARGE	Maximum charging voltage	17V
		Standard charging current	25mA
		Fast charging current	50mA
	DISCHARGE	Range of operating voltage	10V at 17V
		Min tension in discharge	10V
		Max discharge current	125mA
		Lifespan 80% DOD (0,5 C)	500 Cycles
	MAINTENANCE	Frequency of maintenance charges at 20°C	2 Months
	CONTROL ELECTRONICS	Electrical protection	No
		Low voltage power cut	No
		High voltage power cut	No
		Max power cut voltage	No

These devices not only designed to protect the pack in case of an equipment failure. They must not be used to control the discharge. The protection circuits have a response time of a few milliseconds.

MECHANICAL CHARACTERISTICS	Dimensions (+/- 2mm)	Length	51,3mm
		Width	25,75mm
		Depth	34,5mm
	Weight (+/- 5g)	103,5Gr	
	Mechanical protection	Sleeve	
	Wire length (+/- 10mm)	70mm	
	Terminal	JST XHP-2 2,5mm	

CONDITIONS OF USE, STORAGE, AND TRANSPORT	CONDITIONS OF USE	Charging temperature	0 at +45°C
		Discharge temperature	-10 at +45°C
	CONDITIONS OF STORAGE	Storage temperature	-20 at +45°C
		Level of humidity	65,00 %
		Max storage time	2 Years
	TRANSPORT	UN code	Class 9
		ADR/RID classification	Class 9
		IMDG classification	Class 9
		IATA classification	Class 9

INSTRUCTIONS	COMMISSIONING	<ul style="list-style-type: none"> • Check the batteries and the connectors: wires not damaged, battery not swollen, burnt smell, oxidation of the connectors, leak...
	CHARGE	<ul style="list-style-type: none"> • Use an adequate charger • The battery is warmer during the charge: during the first charge, check that the battery's temperature stays in the temperature operating ranges. • In case of an abnormal heating, stop the charge by unplugging the charger within the realms of possibility, remove the battery from the equipment de l'équipement, have the equipment, the battery and the charger checked by a technician.
	CASE OF NON-WATERPROOF BATTERIES	<ul style="list-style-type: none"> • It is normal to observe a release of gas during the charge and use. Do not smoke. Place in suitable premises. • Open batteries need regular maintenance carried out by a qualified technician.
	CASE OF LITHIUM ION BATTERIES	There is a fire hazard with lithium ion batteries in the following cases: overload, short circuit, charge and use outside the voltage and temperature ranges.
	WARNINGS	<ul style="list-style-type: none"> • Read the instructions of your device. • Only use in compatible devices. • Respect the load and storage conditions. • Do not use if the battery is damaged, do not burn, do not pierce, do not dismantle or modify. The protection circuits protect the battery and the equipment: do not deactivate them.

EXPLANATION OF SYMBOLS		Catalogue reference
		Lot number
		manufacturer's address
		To recycle in a suitable salvage and recycling structure
		Read the product sheet and the instruction manual

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