**MGH00811**

Valid from / en cours de validité depuis le / en espera de validación desde el / in corso di validità dal / gültig seit dem / válida a partir de

28/10/2022

<b>Product Sheet (EN)</b>	<b>2</b>
<b>Fiche Produit (FR)</b>	<b>4</b>
<b>Ficha producto (ES)</b>	<b>6</b>
<b>Scheda Prodotto (IT)</b>	<b>8</b>
<b>Eigenschaften (DE)</b>	<b>10</b>
<b>Ficha do Produto (PT)</b>	<b>12</b>

PRODUCT SHEET

valid from 28/10/2022

IDENTIFICATION OF THE BATTERY FOR MEDICAL DEVICE	Type	Rechargeable battery	 
	Commercial designation	Batterie médicale rechargeable POMPE FRESENIUS VIAL 6V 3Ah JST	
	Reference	MGH00811	
	Basic UDI	3660766RechargeablePC	
	EAN	3660766593986	
	Brand	NX	
	Compatible / Original battery	Compatible	
Packaging	Unitary		

RECOMMENDED USAGE

Follow the instructions and recommendations specific to each model, using the technical instructions and document resources from the devices in which the battery is used

Brands	Equipment	Models	PN
Fresenius Vial	Pompe à perfusion	Optima MS	

IDENTIFICATION	GENERAL TECHNICAL CHARACTERISTICS	Chemistry	NiMH
		Type	SC
		IEC designation	5 HRH23/43-1
		Rated voltage	6V
		Nominal capacity	3Ah
		Internal resistance Ω	<30mΩ

The voltage and the actual capacity in use can be affected by several factors, especially the temperature, the discharge current, the pack's history (ex:use, storage), etc

ELECTRICAL CHARACTERISTICS	CHARGE	Maximum charging voltage	8,5V
		Standard charging current (15h)	300mA
		Fast charging current (2,5h)	3000mA
	DISCHARGE	Range of operating voltage	5V to 8,5V
		Min tension in discharge	5V
		Max discharge current	3000mA
		Lifespan 80% DOD (0,5 C)	> 500 Cycles
	CONTROL ELECTRONICS	Electrical protection	Yes
		Low voltage detection threshold	No
		Over voltage detection threshold	No
		Over current detection threshold	Yes

These devices not only designed to protect the pack in case of an equipment failure. They must not be used to control the discharge. The protection circuits have a response time of a few milliseconds.

MECHANICAL CHARACTERISTICS	Dimensions (+/- 2mm)	Length	68,9mm
		Width	43,1mm
		Depth	46,4mm
	Weight (+/- 5g)	284,2Gr.	
	Mechanical protection	Shrink sleeve	
	Wire length (+/- 10mm)	135mm	
	Terminal	JST PHR-5	

CONDITIONS OF USE, STORAGE, AND TRANSPORT	CONDITIONS OF USE	Charging temperature	0 to +45°C	
		Discharge temperature	-30 to +60°C	
	CONDITIONS OF STORAGE	Storage temperature	0 to +40°C	
		Humidity range	15% - 85%	
		Frequency of maintenance charge at 20°C	6 Months	
			Max storage time	2 Years
	TRANSPORT	UN code	3496	
		ADR/RID classification	Not regulated	
		IMDG classification	Not regulated	
		IATA classification	Not regulated	
		Power capacity	18 Wh	

<b>INSTRUCTIONS</b>	<b>COMMISSIONING</b>	<ul style="list-style-type: none"> <li>• Check the batteries and the connectors: wires not damaged, battery not swollen, burnt smell, oxidation of the connectors, leak...</li> <li>• Respect the polarity</li> <li>• Do a full charge with the adequate charger before the first use</li> </ul>
	<b>CHARGE</b>	<ul style="list-style-type: none"> <li>• Use an adequate charger</li> <li>• The battery is warmer during the charge: during the first charge, check that the battery's temperature stays in the temperature operating ranges.</li> <li>• In case of an abnormal heating, stop the charge by unplugging the charger within the realms of possibility, remove the battery from the equipment de l'équipement, have the equipment, the battery and the charger checked by a technician.</li> </ul>
	<b>CASE OF NON-WATERPROOF BATTERIES</b>	<ul style="list-style-type: none"> <li>• It is normal to observe a release of gas during the charge and use. Do not smoke. Place in suitable premises.</li> <li>• Open batteries need regular maintenance carried out by a qualified technician.</li> </ul>
	<b>CASE OF LITHIUM ION BATTERIES</b>	There is a fire hazard with lithium ion batteries in the following cases: overload, short circuit, charge and use outside the voltage and temperature ranges.
	<b>WARNINGS</b>	<ul style="list-style-type: none"> <li>• Read the instructions of your device.</li> <li>• Only use in compatible devices.</li> <li>• Respect the load and storage conditions.</li> <li>• Do not use if the battery is damaged, do not burn, do not pierce, do not dismantle or modify. The protection circuits protect the battery and the equipment: do not deactivate them.</li> </ul>

<b>EXPLANATION OF SYMBOLS</b>	<b>LABELLING</b>	
		Catalogue reference
		Lot number
		Manufacturer's address
		To recycle in a suitable salvage and recycling structure
		Read the product sheet and the instruction manual
		Production date
		Fragile, handle with care . Can be broken or damaged if not handled carefully
		Keep away from sunlight. "Keep away from heat".
		Keep dry. Needs to be protected from moisture. "Keep away from rain"
		Temperature limit. Indicates the temperature limits to which the medical device can be safely exposed. The upper and lower limits of temperature are indicated adjacent to the upper and lower horizontal lines.
		Humidity limitation. Indicates the range of humidity to which the medical device can be safely exposed. The humidity limitation are indicated adjacent to the upper and lower horizontal lines.
		Unique device identifier. Indicates a carrier that contains unique device identifier information NOTE This symbol identifies the UDI carrier, including the AIDC (Automatic Identification and Data Collection) and human readable information.

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