MGL00659

Valid from / en cours de validité depuis le / en espera de validación desde el / in corso di validità dal / gültig seit dem / válida a partir de

02/08/2023

Product Sheet (EN)	2
Fiche Produit (FR)	4
Ficha producto (ES)	6
Scheda Prodotto (IT)	8
Eigenschaften (DE)	10
Ficha do Produto (PT)	12

PRODUCT SHEET

valid from 02/08/2023

IDENTIFICATION OF THE BATTERY FOR MEDICAL DEVICE	Type	Rechargeable battery	 
	Commercial designation	Batterie médicale rechargeable Invacare Aquatec & Orca 11.1V 2.2Ah COSSE	
	Reference	MGL00659	
	Basic UDI	3660766RechargeablePC	
	EAN	3660766621788	
	Brand	NX, a brand of Upergy	
	Compatible / Original battery	Compatible	
Packaging	Unitary		

RECOMMENDED USAGE

Follow the instructions and recommendations specific to each model, using the technical instructions and document resources from the devices in which the battery is used

Brands	Equipment	Models	PN
Invacare	Bath lift	Aquatec Orca	TE-1561-1

IDENTIFICATION	GENERAL TECHNICAL CHARACTERISTICS	Chemistry	Li-Ion
		Type	18650
		IEC designation	3 ICR19/66-1
		Rated voltage	11,1V
		Nominal capacity	2,2Ah
	Internal resistance Ω	<105mΩ	

The voltage and the actual capacity in use can be affected by several factors, especially the temperature, the discharge current, the pack's history (ex:use, storage), etc

ELECTRICAL CHARACTERISTICS	CHARGE	Maximum charging voltage	12,6V
		Standard charging current (15h)	220mA
		Fast charging current (2,5h)	2200mA
	DISCHARGE	Range of operating voltage	7,5V to 12,6V
		Min tension in discharge	7,5V
		Max discharge current	2200mA
		Lifespan 80% DOD (0,5 C)	> 500 Cycles
	CONTROL ELECTRONICS	Electrical protection	Yes
		Low voltage detection threshold	No
		Over voltage detection threshold	No
	Over current detection threshold	Yes	

These devices not only designed to protect the pack in case of an equipment failure. They must not be used to control the discharge. The protection circuits have a response time of a few milliseconds.

MECHANICAL CHARACTERISTICS	Dimensions (+/- 2mm)	Length	36,6mm
		Width	34,2mm
		Depth	66mm
	Weight (+/- 5g)	139g	
	Mechanical protection	None	
	Wire length (+/- 10mm)	70mm	
	Terminal	Cosses	

CONDITIONS OF USE, STORAGE, AND TRANSPORT	CONDITIONS OF USE	Charging temperature	0 to +40°C
		Discharge temperature	-20 to +70°C
	CONDITIONS OF STORAGE	Storage temperature	0 to +30°C
		Humidity range	30% - 75%
		Frequency of maintenance charge at 20°C	6 Months
		Max storage time	2 Years
	TRANSPORT	UN code	3480
		ADR/RID classification	Exempt
		IMDG classification	Exempt
		IATA classification	Dangerous
	Power capacity	24,42 Wh	

INSTRUCTIONS	COMMISSIONING	<ul style="list-style-type: none"> • Check the batteries and the connectors: wires not damaged, battery not swollen, burnt smell, oxidation of the connectors, leak... • Respect the polarity • Do a full charge with the adequate charger before the first use
	CHARGE	<ul style="list-style-type: none"> • Use an adequate charger • The battery is warmer during the charge: during the first charge, check that the battery's temperature stays in the temperature operating ranges. • In case of an abnormal heating, stop the charge by unplugging the charger within the realms of possibility, remove the battery from the equipment de l'équipement, have the equipment, the battery and the charger checked by a technician.
	CASE OF NON-WATERPROOF BATTERIES	<ul style="list-style-type: none"> • It is normal to observe a release of gas during the charge and use. Do not smoke. Place in suitable premises. • Open batteries need regular maintenance carried out by a qualified technician.
	CASE OF LITHIUM ION BATTERIES	There is a fire hazard with lithium ion batteries in the following cases: overload, short circuit, charge and use outside the voltage and temperature ranges.
	WARNINGS	<ul style="list-style-type: none"> • Read the instructions of your device. • Only use in compatible devices. • Respect the load and storage conditions. • Do not use if the battery is damaged, do not burn, do not pierce, do not dismantle or modify. The protection circuits protect the battery and the equipment: do not deactivate them.

EXPLANATION OF SYMBOLS	LABELLING	
		Catalogue reference
		Lot number
		Manufacturer's address
		To recycle in a suitable salvage and recycling structure
		Read the product sheet and the instruction manual
		Production date
		Fragile, handle with care . Can be broken or damaged if not handled carefully
		Keep away from sunlight. "Keep away from heat".
		Keep dry. Needs to be protected from moisture. "Keep away from rain"
		Temperature limit. Indicates the temperature limits to which the medical device can be safely exposed. The upper and lower limits of temperature are indicated adjacent to the upper and lower horizontal lines.
		Humidity limitation. Indicates the range of humidity to which the medical device can be safely exposed. The humidity limitation are indicated adjacent to the upper and lower horizontal lines.
		Unique device identifier. Indicates a carrier that contains unique device identifier information NOTE This symbol identifies the UDI carrier, including the AIDC (Automatic Identification and Data Collection) and human readable information.

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