



MEDICAL BATTERY 11.1V 6.8Ah

BATTERIE MÉDICALE SONOSITE POUR S-NERVE 10.8V 5.8Ah

MAIN INFORMATION / INFORMATIONS GÉNÉRALES

BATTERY TYPE	TYPE DE BATTERIE	Compatible battery/ Batterie compatible
BRAND	MARQUE	Mindray
TECHNOLOGY	TECHNOLOGIE	LI-ION
VOLTAGE	TENSION	11.1V
CAPACITY	CAPACITÉ	6.8Ah

ADDITIONAL INFORMATION / INFORMATIONS COMPLÉMENTAIRES

UNIT PACKING	CONDITIONNEMENT UNITAIRE	
---------------------	--------------------------	--

COMPATIBLE PRODUCTS / PRODUITS COMPATIBLES*

PART NUMBER / PART NUMBER	0146-00-0069
	0146-00-0099
	115-018011-00
	115-018015-00
	11501801500
	6132
	9201-30-35944
	EE02100
	LI23S001A
	M05-0100001-06
	M05010000106
	MDR-M05-010001-06
MANUFACTURER – MODEL / MARQUES – MODÈLES	
DATASCOPE	ACCUTORR PLUS
DATASCOPE	ACCUTORR V
DATASCOPE	DATASCOPE EX 0146-00-0099 - TRIO
DATASCOPE	DPM
DATASCOPE	DPM3
DATASCOPE	DPM4
DATASCOPE	DPM5
DATASCOPE	MONITEUR ACCUTOR +
DATASCOPE	MONITEUR DE TRANSPORT PASSPORT 2
DATASCOPE	MONITEUR SPECTRUM
DATASCOPE	MONITEUR TRIO
DATASCOPE	PASSPORT PM7000
DATASCOPE	PASSPORT PM8000
DATASCOPE	PASSPORT PM9000
DATASCOPE	PASSPORT V
DATASCOPE	PASSPORT 3
DATASCOPE	PM-8000 EXPRESS
DATASCOPE	PM-9000
DATASCOPE	PM-9000 EXPRESS
DATASCOPE	SPECTRUM OR
DATASCOPE	VS-800
DATASCOPE	VS800



* Non exhaustive list / Liste non-exhaustive

www.enix-energies.com • Date: 2021-06-25



OUTLINE SAFETY WARNING: USE ONLY WITH IN THE ALLOWED PARAMETERS. Do not short circuit or over-load the battery. Charge only using an approved charger designed specifically to charge this battery. Do not heat above maximum temperatures indicated. Never crush, mutilate, puncture or abuse the battery. Do not dismantle the pack or disable any of the protective devices or circuits. DO NOT USE THE BATTERY IF YOU SUSPECT IT MAY BE FAULTY OR DAMAGED.

© Copyright Enix Energies 2005. NB: This document and the product design are the intellectual property of Enix Energies. No document or design may be copied or used for commercial purposes without written permission of Enix Energies. Users must satisfy themselves, by means of testing etc, that products are suitable for their application. Data in this document is for general guidance only; consult cell manufacturers data for definitive information. Information is given free of charge and in good faith, but no responsibility can be accepted for any errors or omission or costs or losses or liabilities arising from the use of this information. All business is conducted to Enix Energies terms and conditions only.